

## Injection Areas

### dermaren™ FINE with Lidocaine


- Fine Lines
- Crow's Feet
- Glabella Lines
- Forehead Lines
- Earlobe
- Perioral Wrinkles

### dermaren™ DEEP with Lidocaine

- ★ Deep Lines/ Augmentation
- ★ Nasolabial Lines
- ★ Glabella Lines
- ★ Forehead Lines
- ★ Laugh Lines
- ★ Cheek/ Chin/ Forehead Augmentation
- ★ Lip Contouring
- ★ Nasal Augmentation

### dermaren™ SUB-Q with Lidocaine

- Deep Lines/ Augmentation
- Nasolabial Lines
- Laugh Lines
- Cheek/ Chin/ Forehead Augmentation
- Lip Contouring
- Nasal Augmentation
- Breast Augmentation



# dermaren™

CROSS LINKED HYALURONIC ACID  
With Lidocaine

*Long Lasting Confidence*



100% Cross-linked Hyaluronic Acid Filler

*Acmedi* KOREA

Zip. 07997, #611, Hyundai 41 Tower, 293, Mokdongdong-ro, Yangcheon-gu, Seoul, Korea  
T. +82 70 7590 3883 F. +82 2 6971 9764 E. sales@acmedikorea.com  
www.acmedikorea.com

*Acmedi* KOREA

## About dermaren™

**dermaren™** is a Non-animal Cross-linked Hyaluronic Acid Soft Tissue Filler.

**dermaren™** has proved its Effectiveness and Safety through several clinical trials compared with Restylane® (Qmed, Sweden). 1)

**dermaren™** keeps skin moisturized and stimulates the aged cells to generate.

**dermaren™** is an ideal product for volumizing effect and high patients' satisfaction.

## Various Choices of Product Types


	 dermaren™ FINE with Lidocaine	 dermaren™ DEEP with Lidocaine	 dermaren™ SUB-Q with Lidocaine
			
Composition	HA 24mg/ml Lidocaine 0.3%	HA 24mg/ml Lidocaine 0.3%	HA 24mg/ml Lidocaine 0.3%
Suggested Uses	Superficial Lines : Worry lines, Crow's feet, Neck wrinkles, Periorbital lines, Perioral lines	Mid/ Deep Folds and Wrinkles : Nasolabial folds, Cheeks, Chins, Lips,	Deep Folds and Wrinkles : Facial contours, Cheeks, Lips, Nose enhancement
Injection Level	Superficial dermis	Mid/ Deep dermis	Deep dermis Subcutaneous
Syringe Volume	1.1ml	1.1ml	1.1ml
Needle Size	30G	27G	25G
Storage	2-25 °C		
Duration	12-18months		
BDDE	Not Detected		
Shelf Life	24months		

\*1) Clinical Study Results of Across Co., Ltd

## Features 1) Monophasic Blob Structure


**dermaren™** has Regular and Dense Lattice comparing to other products which allows Consistent and Smooth injection throughout the treatment and helps naturally harmonized look.

dermaren™




Product [ J ]


Product [ P ]



## Features 2) Extrusion Force & Injection Pain Comparison



**dermaren™**, its consistent extrusion force enables doctors to control the injections easily and safely throughout the procedure.



**dermaren™** contains 0.3% Lidocaine, so it minimizes the pain during and or post injection for the patients.