



PRODUCT CATALOG

UCARETRON INC.






Wireless Massager



Who use it?



Workers



Housewives



Students



Drivers



Seniors

EASY CONTROL

The massager is connected to a smart phone by Bluetooth and controlled by app(Android/iOS)

As the massager is often put on the shoulder or back, wireless control is inevitable



User can operate two massagers simultaneously

They can get a massage on different parts of their body at the same time just like hand-massage



Each massager is separately controlled by time-delay, stimulation level.

Music Beat Massage



The massager stimulates to the music beat.
As our technology changes the beat to the waveform of
massage,
the user can take a massage more dynamically.



40 songs are pre-saved in the app

- Intensive beat music : 20 songs
- Piano : 10 song
- Meditation music : 10 songs


Massage Mode




Stimulation Level – Up to 20



Various Pad



Various Pad




HYDROGEL ON PAD

Hydrogel on the pad
maximizes effectiveness of
electrical stimulation.





If the adhesion is weak, it can be used
more after being washed in the water

Battery



PRODUCT COMPOSITION – BEATCONNECT BASIC PACKAGE



Large Pad x 1pcs



Asymmetric Pad x 1pcs





Small Pad x 1pcs



Wrist Pad x 1pcs

PRODUCT COMPOSITION – BEATCONNECT STEREO PACKAGE



Large Pad x 2pcs



Asymmetric Pad x 2pcs




Small Pad x 2pcs




Wrist Pad x 2pcs

FOOT PAD – OPTIONAL ACCESSORY



SPECIFICATION

Model No	UM-200B	
Input	5V	
Current Consumption	20mA	
Battery	Li-ion Polymer, 3.7V / 70mAh	
Material	ABS / Polyurethane(Pad)	
Size(mm)	Beat Connected: 40x40 <ul style="list-style-type: none">• Small Pad : 146x57• Large Pad : 195x68• Asymmetric Pad : 177x70,• Wrist Pad : 175x35	
Weight	Beat Connected : 15.2g <ul style="list-style-type: none">• Small Pad : 14g• Large Pad : 27.6g• Asymmetric Pad : 22.6g• Wrist Pad : 14.2g	

Check baby's body temperature all the times



T-Patch

Color-indicating Thermometer



Sticker Thermometer

- Simply attached on the skin
- Monitoring color change
- Checking body temperature



PACKAGE

Medical device
의료기기

Body temperature is indicated by color change of sticker.
스티커의 색상으로 체온을 표시합니다.



< 35°C 35°C ~ 36.5°C > 37.5°C

**티패치
칼라 스티커 온도계**

T-Patch Sticker Thermometer

- Simple attachment to skin
피부에 간단히 부착
- Monitoring color change
색상 변화 관찰
- Measuring body temperature anytime
언제든지 체온 관찰 가능

20ea (20개입)

U (주)유케어드론
PATENT

Patent pending (특허 10-2020-0127739)
Utility pending (실용신안 20-2020-0003872)

20 Stickers

